

QUALITY CONTROL CERTIFICATE

Cat. No. MAB602P	Monoclonal Antibody to Respiratory Syncytial Virus
-------------------------	---

Clone Determination	12D11
Immunogen	Respiratory Syncytial Virus
Immunoglobulin Class	IgG _{2a}
Specificity	Reacts with RSV protein of 33kDa and ~190kDa as determined by Western Blot
Purification	Protein G Purification from Supernatant
Lot Number	
Concentration	1mg/mL
Buffer	0.15M PBS
PH	7.0-7.5
Preservatives	None
Expiration date	
Quantity	
Volume	
Applications	ELISA: Cell Substrate IFA Western Blot ELISA as a Capture or Tracer Chromatographic Immunoassay
Working dilution	Determined by User
Storage	Store at -20C to -80C. Avoid repeat freeze- thaw cycles

FOR RESEARCH AND *IN VITRO* USE ONLY

MAINE BIOTECHNOLOGY SERVICES, INC.
1037R FOREST AVENUE PORTLAND, ME 04103 USA
Tel.: 800-925-9476 / 207-797-5454 / Fax: 207- 797-5595