

QUALITY CONTROL CERTIFICATE

Cat. No. MAB838P	Monoclonal Antibody to <i>Giardia lamblia</i>
-------------------------	--

Clone Determination	6261
Immunogen	Purified GSA65 antigen
Immunoglobulin Class	IgG2b
Specificity	Specific to the 65 kDa <i>Giardia</i> antigen, detects cysts.
Purification	Purified by Protein A chromatography from ascites.
Lot Number	
% Purity	> 90%
Concentration	
Buffer	0.01M PBS, pH 7.2
Preservatives	0.1% NaN ₃
Expiration date	< 6 months at 4C or 2 years in aliquots at -20C
Quantity	
Volume	
Applications	ELISA, IFA, WB
Working dilution	Determined by user
Storage	< 6 months at 4C or 2 years in aliquots at -20C

FOR RESEARCH AND *IN VITRO* USE ONLY

MAINE BIOTECHNOLOGY SERVICES, INC.
1037R FOREST AVENUE PORTLAND, ME 04103 USA
Tel.: 800-925-9476 / 207-797-5454 / Fax: 207- 797-5595