

**QUALITY CONTROL CERTIFICATE**

<b>Cat. No. MAB839P</b>	<b>Monoclonal Antibody to <i>Giardia lamblia</i></b>
-------------------------	------------------------------------------------------

<b>Clone Determination</b>	6271
<b>Immunogen</b>	Purified <i>Giardia lamblia</i> cysts
<b>Immunoglobulin Class</b>	IgG3
<b>Specificity</b>	Reactive with <i>Giardia</i> cysts in IFA and ELISA. Non-reactive with <i>Cryptosporidium</i> and <i>Microsporidia</i> cysts
<b>Purification</b>	Purified by Protein A chromatography from ascites.
<b>Lot Number</b>	
<b>% Purity</b>	> 90%
<b>Concentration</b>	
<b>Buffer</b>	0.01M PBS, pH 7.2
<b>Preservatives</b>	0.1% NaN <sub>3</sub>
<b>Expiration date</b>	< 6 months at 4C or 2 years in aliquots at -20C
<b>Quantity</b>	
<b>Volume</b>	
<b>Applications</b>	ELISA, IFA
<b>Working dilution</b>	Determined by user
<b>Storage</b>	< 6 months at 4C or 2 years in aliquots at -20C

**FOR RESEARCH AND *IN VITRO* USE ONLY**

**MAINE BIOTECHNOLOGY SERVICES, INC.**  
1037R FOREST AVENUE PORTLAND, ME 04103 USA  
Tel.: 800-925-9476 / 207-797-5454 / Fax: 207- 797-5595

email:sales@mainebiotechnology.com www.mainebiotechnology.com

rev. 12/01/08