

**QUALITY CONTROL CERTIFICATE**

<b>Cat. No. MAB840P</b>	<b>Monoclonal Antibody to <i>Campylobacter jejuni</i></b>
-------------------------	---

<b>Clone Determination</b>	7711
<b>Immunogen</b>	Native <i>C. jejuni</i> organism
<b>Immunoglobulin Class</b>	IgG1
<b>Specificity</b>	Appears specific to <i>C. jejuni</i> via whole cell ELISA. Reactive with an ~100 kDa antigen in WB. Non-reactive with <i>H. pylori</i> , <i>E. coli</i> , <i>Salmonellae</i> , <i>Shigellae</i> , <i>P. aeruginosa</i> , <i>Yersinia</i> and <i>Citrobacter</i> .
<b>Purification</b>	Purified by Protein A chromatography from ascites.
<b>Lot Number</b>	
<b>% Purity</b>	> 90%
<b>Concentration</b>	
<b>Buffer</b>	0.01M PBS, pH 7.2
<b>Preservatives</b>	0.1% NaN <sub>3</sub>
<b>Expiration date</b>	< 6 months at 4C or 2 years in aliquots at -20C
<b>Quantity</b>	
<b>Volume</b>	
<b>Applications</b>	ELISA, WB
<b>Working dilution</b>	Determined by user
<b>Storage</b>	< 6 months at 4C or 2 years in aliquots at -20C

**FOR RESEARCH AND *IN VITRO* USE ONLY**

**MAINE BIOTECHNOLOGY SERVICES, INC.**  
1037R FOREST AVENUE PORTLAND, ME 04103 USA  
Tel.: 800-925-9476 / 207-797-5454 / Fax: 207- 797-5595