

QUALITY CONTROL CERTIFICATE

Cat. No. MAB841P **Monoclonal Antibody to Campylobacter jejuni**

Clone Determination	7721
Immunogen	Native <i>C. jejuni</i> organism
Immunoglobulin Class	IgG1
Specificity	Appears specific to <i>C. jejuni</i> via whole cell ELISA. Reactive with a ~100 kDa antigen in WB. Non-reactive with <i>H. pylori</i> , <i>E. coli</i> , <i>Salmonellae</i> , <i>Shigellae</i> , <i>P. aeruginosa</i> , <i>Yersinia</i> and <i>Citrobacter</i> .
Purification	Purified by Protein A chromatography from ascites.
Lot Number	
% Purity	> 90%
Concentration	
Buffer	0.01M PBS, pH 7.2
Preservatives	0.1% NaN ₃
Expiration date	< 6 months at 4C or 2 years in aliquots at -20C
Quantity	
Volume	
Applications	ELISA, WB
Working dilution	Determined by user
Storage	< 6 months at 4C or 2 years in aliquots at -20C

FOR RESEARCH AND *IN VITRO* USE ONLY

MAINE BIOTECHNOLOGY SERVICES, INC.
1037R FOREST AVENUE PORTLAND, ME 04103 USA
Tel.: 800-925-9476 / 207-797-5454 / Fax: 207- 797-5595