

### **Which animal should I choose?**

The following questions can help you make an optimum selection of species for your project:

- How much antibody do I need for the project development and final application?
- Does the host animal make a difference to the final application?
- Will the antibody need to be affinity purified or Ig purified?
- How much antigen do I have available?

### **Expected Antibody Yield Guideline:**

**Rabbits:** Neat serum antibody yield is generally 4-10 mg/ml of immunoglobulin.  
Specific antibody yield is generally 1-2 mg/ml.

**Goats & Sheep:** Neat serum antibody yield is generally 5-15 mg/ml of immunoglobulin.  
Specific antibody yield is generally 2-3 mg/ml.